

# Infectious Disease Portfolio: Global Challenges & Illumina Solutions

---

Tung, Vu

Senior FAS, Vietnam

09 June 2024

---

OUR MISSION

Improve human health  
by unlocking the power of  
the genome

---



# Whom we serve

## KEY MARKETS



Reproductive  
health



Genetic  
disease



Oncology



Research /  
Applied genomics



Infectious  
disease



Drug  
discovery



Multomics

## CUSTOMERS



Universities  
and academic  
research centers



Pharmaceutical  
companies



Genome  
centers



Biotechnology  
companies



Hospitals



Consumer  
genetics  
companies



Government  
agencies

# The sequencing workflow

## 1. Library preparation



A sequencing library contains the DNA of interest to the experiment.

Adapters are added to both ends of the molecules to uniquely identify them.

Many libraries may be pooled together for a single sequencing run.

## 2. Cluster growth/generation



Libraries are loaded into a flow cell, a glass slide where sequencing chemistry occurs.

Each DNA fragment is amplified into distinct clonal clusters through bridge amplification.

## 3. Sequencing



Fluorescently labeled nucleotides are incorporated with complementary DNA bases of interest.

The instrument excites the fluorescent labels and captures an image of the flow cell, allowing identification.

## 4. Data analysis



Primary analysis, on the instrument, converts the images into base calls.

Secondary analysis generates alignments and variant detection.

Tertiary analysis allows for annotation, filtering, and interpretation.

# More than 21,000 systems installed around the world\*

## Our sequencing and array instruments

### High-throughput



NovaSeq™ X  
NovaSeq™ X Plus



NovaSeq™ 6000



iScan™ System



NextSeq™ 1000  
NextSeq™ 2000



NextSeq™ 550



NextSeq™ 500

### Low-throughput



iSeq™ 100



MiniSeq™



MiSeq™



NovaSeq™ 6000Dx†



NextSeq™ 550Dx†



MiSeq™ Dx†

### IVD instruments

\*Net installed base as of fiscal year 2023

# Sequencing systems and key applications overview

		Large whole-genome sequencing (WGS) <i>(human, plant, animal)</i>	Small WGS <i>(microbe, virus)</i>	Exome sequencing	Targeted gene sequencing	Whole-transcriptome sequencing	Gene expression profiling with mRNA-Seq	Targeted gene expression profiling
High-throughput	<b>NovaSeq X / NovaSeq X Plus</b> <ul style="list-style-type: none"> <li>Our most powerful sequencer, equipped with XLEAP-SBS chemistry</li> </ul>							
	<b>NovaSeq 6000</b> <ul style="list-style-type: none"> <li>Scalable throughput for broad, deep sequencing</li> </ul>							
Mid-throughput	<b>NextSeq 1000/2000</b> <ul style="list-style-type: none"> <li>Flexibility, affordability, and scalability</li> </ul>							
	<b>NextSeq 500/550</b> <ul style="list-style-type: none"> <li>Tunable output with sequencing and array capabilities</li> </ul>							
Low-throughput	<b>MiSeq</b> <ul style="list-style-type: none"> <li>Our first benchtop sequencer</li> </ul>							
	<b>MiniSeq</b> <ul style="list-style-type: none"> <li>&lt; 1 day turnaround time</li> </ul>							
	<b>iSeq 100</b> <ul style="list-style-type: none"> <li>Our most affordable sequencer</li> </ul>							

Note: Only key applications are highlighted, which does not reflect each system's entire set of capabilities.

**For Research Use Only. Not for use in diagnostic procedures.**

- 
- 1 Introduction**
  - 2 Portfolio &. Methods**
  - 3 Single Pathogen**
  - 4 Multi-pathogen & Syndromic Panels**
  - 5 Discovery**



# Global Challenges in Infectious Disease



## Outbreak, Pandemic, Epidemic Disease Management

Many new factors contribute to an increase in the transmissibility and severity of infectious diseases



## Spread of Antimicrobial Resistance 'AMR'

Increased threat due to multi-drug resistant organisms through the environment



## Diagnostic for Patient Management

Access to low cost, accurate IVD's is still a major challenge in many under-developed countries



## Infectious Diseases by the Numbers

770M	Cumulative cases of <b>COVID-19</b> <sup>1</sup>
400K	Respiratory deaths each year caused by <b>influenza virus</b> <sup>2</sup>
10.6M	People that fell ill from <b>tuberculosis</b> in 2021 <sup>3</sup>
39M	People living with <b>HIV</b> in 2022 <sup>4</sup>
3.8%	Global GDP lost due to <b>antimicrobial resistance</b> by 2050 <sup>5</sup>
247M	<b>Malaria</b> cases in 84 malaria-endemic countries in 2021 <sup>3</sup>
1.5M	New <b>hepatitis B</b> infections globally in 2019 <sup>3</sup>
\$28-45B	Direct annual cost of treating <b>healthcare associated infections</b> in the US <sup>6</sup>

1. <https://covid19.who.int/> (Accessed Sept. 7, 2023); 2. <https://ourworldindata.org/influenza>; 3. World health statistics 2023: monitoring health for the SDGs, Sustainable Development Goals. Geneva: World Health Organization; 2023.; 4. <https://www.unaids.org/en/resources/fact-sheet>; 5. Global antimicrobial resistance and use surveillance system (GLASS) report 2022. Geneva: World Health Organization; 2022.; 6. Gidey K, et al. Clinical and economic burden of healthcare-associated infections: A prospective cohort study. PLoS One. 2023 Feb 23;18(2):e0282141.



Infectious  
diseases impact  
human health  
and economics

# Illumina technologies – a unique advantage against infectious diseases




Flexible and scalable solutions to enable accessibility worldwide

- ✓ Novel pathogen ID and drug resistance characterization
- ✓ Insights for clinical assay development
- ✓ Pathogen monitoring and genomic surveillance
- ✓ Host genetics and immune response
- ✓ Drug and vaccine development

# Portfolio & Methods


## Introduction

---




# Enabling the ID&M Portfolio with End-to-End Workflows



Library Prep/Panel	Systems	Current Software- on BSSH
Illumina DNA Prep		DRAGEN™ Metagenomics & SRST2
COVIDSeq		DRAGEN™ COVID Lineage
Illumina and GenoScreen Deeplex Myc-TB Combo Kit		Deeplex Web App
Illumina Microbial Amplicon Prep (IMAP) Illumina Microbial Amplicon Prep – Influenza		DRAGEN™ Targeted Microbial
Respiratory Virus Enrichment Kit (RVEK, formerly RVOP)		DRAGEN™ RNA Pathogen Detection
Viral Surveillance Panel (VSP) Pan-Coronavirus Panel		DRAGEN™ Microbial Enrichment
Respiratory Pathogen ID/AMR Panel (RPIP) Urinary Pathogen ID/AMR Panel (UPIP)		Explify™
Ribo-Zero Plus Microbiome		Microbiome Metatranscriptomics

# NGS Methods Enable Genomic Investigation



# Single Pathogen

“I know what I am looking for...”

---





# Illumina DNA Prep

Illumina library prep with the flexibility for whole-genome sequencing, amplicon & shotgun metagenomics methods for nucleotide resolution.



## Applications

Genomic surveillance (HAI, Foodborne pathogens), gene-specific amplicons (16S, AMR genes), pathogen discovery, microbiome analysis & host/pathogen studies and more

## Methods

Whole genome sequencing, amplicon & shotgun metagenomics or metatranscriptomics

## Other technologies

Microbial culture, Sanger sequencing, qPCR (presence/absence), pulse field electrophoresis (PFGE)

## Features

Fastest Illumina library prep workflow, with ~3.5 hours total time

Flexibility to accommodate variations in sample types (gDNA, blood, saliva, dried blood spots, and direct from culture), DNA input amount (1-500 ng), and methods (WGS, amplicon or metagenomics)

Optimized library prep performance; delivers consistent insert sizes, uniform coverage, and optimized performance, regardless of DNA input amount or genome size.

The bead-based technology minimizes bias and opportunities for error, resulting in highly reproducible sequencing data.

Compatible with all Illumina platforms

Automatable through our [many partners](#)

## Current users:

- Researchers, Academic labs
- Public health labs
- Clinical research
- Applied Microbiology Labs (Pharma, manufacturing, direct-to-consumer)
- Service providers

# Amplicon Targeted Assays



01

Amplicon is faster and easier than hybrid capture

02

Amplicon-based COVIDSeq™ familiar to many customers

03

Risk of mutations causing low coverage, primer design may need updating over time

GenoScreen Deeplex Myc-TB	Identification of MTBC, 15 anti-TB drugs, and >100 non-TB mycobacteria
COVIDSeq	Full genome SARS-CoV-2
Illumina Microbial Amplicon Prep – Influenza (IMAP-Flu)	Full genome Influenza A&B
Illumina Microbial Amplicon Prep (IMAP)	Open platform for ID&M research & surveillance applications

# Illumina and GenoScreen Deeplex Myc-TB Combo Kit

Culture-free amplification of *Mycobacterium tuberculosis* drug targets with strain identification, as well as Nontuberculous Mycobacteria species identification

\*Detection of heteroresistance down to 1-3% will depend on coverage depth



For Research Use Only. Not for use in diagnostic procedures.

## Applications

- Tuberculosis genomic surveillance
- Anti-tuberculosis drug resistance
- Mycobacterium tuberculosis characterization
- Non-Mycobacterium tuberculosis identification

## Methods

Amplicon Sequencing (tNGS)

## Other technologies

Microbial culture, qPCR (presence/absence), Line Probe Assay

## Features

Culture-free or cultured isolates

Prediction of resistance to 15 anti-TB drugs in 24-48h

Detection of minority variants potentially conferring heteroresistance down to 1-3%

MTBC strain typing (spoligotyping and SNP typing)

Identification of more than 100 Nontuberculous Mycobacteria (NTM)

Secure and automated analysis and easy interpretation of results via Deeplex Web App

Compatible with iSeq™ 100, MiniSeq™, MiSeq™, and NextSeq™ 550 platforms

## Current users:

- National TB Programs, Public health or Ministry of Health labs
- Clinical research
- Applied Microbiology Labs (Pharma, manufacturing, direct to consumer)
- Researchers, Academic labs
- Service providers

# WHO guidance of products that meet tNGS class-based performance criteria

3 tNGS solutions that met the inclusion criteria for WHO evaluation

Where a product has not yet met the requirements for a specific drug (i.e. no check mark below), further improvements to the product and WHO review of the evidence will be necessary before clinical use.

Drug class	TB Drug	GenoSreen Deeplex Myc-TB Combo kit	Company O	Company S	Company A
First-line	rifampicin	✓	✓		<b>Company A Assay</b> did not meet the inclusion criteria and as such, is currently not WHO endorsed until further data can be reviewed by WHO
	isoniazid	✓	✓		
	fluoroquinolones	✓	✓		
	ethambutol	✓		✓	
	pyrazinamide	✓			
Second-line	bedaquiline	✓			
	linezolid	✓	✓		
	clofazimine	✓			
	amikacin	✓	✓		
	streptomycin	✓	✓		

WHO Rapid Communication: Use of targeted next-generation sequencing to detect drug-resistant tuberculosis. July 2023 <https://www.who.int/publications/i/item/9789240076372>

# COVIDSeq Assay and COVIDSeq Test

Whole genome amplicon design for the detection  
and characterization of the SAR-CoV-2 RNA virus



For Research Use Only. Not for use in diagnostic procedures.

## Applications

- SARS-CoV-2 identification
- Novel mutations
- SARS-CoV-2 Surveillance
- Vaccine efficacy

## Methods

Amplicon Sequencing (WGS)

## Other technologies

Sanger sequencing, qPCR (presence/absence), lateral flow assays, shotgun metagenomics

## Features

Full genome design of SARS-CoV-2 with 99 targets including S-genes, for highly accurate detection.

Flexible input for both nasal swabs and wastewater samples

Analysis powered by DRAGEN™ COVID Lineage App

Compatible with all Illumina Systems

COVIDSeq Assay (96 samples)

COVIDSeq Test (3072 samples)

## Current user:

- Public health labs
- Researchers, Academic labs
- Clinical research
- Applied Microbiology Labs (Pharma, manufacturing, direct to consumer)
- Service providers

\*not tested on NovaSeqX

# Illumina Microbial Amplicon Prep Influenza (IMAP-Flu)

Whole genome amplicon design for the detection and characterization of Influenza A/B RNA virus



For Research Use Only. Not for use in diagnostic procedures.

## Applications

- Genomic surveillance
- Influenza virus identification
- Seasonal virus characterization

## Methods

Amplicon Sequencing

## Other technologies

Sanger sequencing, qPCR (presence/absence), PCR.

## Features

Amplicon-based library preparation kit with reagents and indexes for 48 samples with data analysis included

Flexibility to support multiple sample types (including nasal swabs and isolates)

>98% coverage for both influenza A and B viruses

Works with multiple subtypes including H1N1, H3N2, and H5N1 (highly pathogenic avian influenza)

Streamlined workflow, 9 hours from extracted RNA to library

Analysis powered by DRAGEN Targeted Microbial App

Compatible with all Illumina Systems\*, ideal for benchtop systems

## Current users:

- Researchers, Academic labs
- Public health labs
- Clinical research
- Applied Microbiology Labs (Pharma, manufacturing, direct-to-consumer)
- Service providers

\*not tested on NovaSeqX

# Illumina Microbial Amplicon Prep (IMAP)

Unique and flexible amplicon library prep solution allowing researchers to use lab designed or commercially available primers for targeted sequencing of viruses and select markers from larger microbial species



For Research Use Only. Not for use in diagnostic procedures.

## Applications

- Genomic surveillance
- Viral surveillance
- Functional Genomics
- Recombinant virus sequencing
- Vaccine efficacy and development
- Wastewater surveillance

## Methods

Amplicon Sequencing

## Other technologies

Amplicon Sequencing  
Sanger sequencing, qPCR (presence/absence), PCR.

## Features

Amplicon-based library preparation kit with reagents and indexes for 48 samples with data analysis included

Based on COVIDSeq Assay chemistry and protocol

Use with lab-designed primers, published primers, or commercially available primer sets

Starting from extracted RNA, DNA, or total nucleic acid inputs

Flexibility to accommodate variations in sample types (e.g. swabs, wastewater, culture)

Analysis powered by DRAGEN Targeted Microbial App

Compatible with all Illumina Systems\*, ideal for benchtop systems

## Current users:

- Researchers, Academic labs
- Public health labs
- Clinical research
- Applied Microbiology Labs (Pharma, manufacturing, direct to consumer)
- Service providers

\*not tested on NovaSeqX



# Illumina Microbial Amplicon Prep

## Tested & Demonstrated Protocols

\* The list is for informational purposes only and does not imply the protocols are validated or supported by Illumina.



For Research Use Only. Not for use in diagnostic procedures.

Illumina Tested <sup>1</sup>	Customer Demonstrated
Chikungunya	<a href="#">Dengue I-IV</a> <sup>2</sup>
RSV	<a href="#">Lumpy skin disease</a> <sup>3</sup>
SARS-CoV-2: ARTIC v4	<a href="#">Mpox</a> <sup>4</sup>
Zika	<a href="#">RSV</a> <sup>5</sup>
	SARS-CoV-2: ARTIC v4.1 <sup>6</sup> & v5.3.2 <sup>7</sup>

1. Data on file
2. <https://dx.doi.org/10.17504/protocols.io.kqdg39xxeg25/v2>
3. <https://doi.org/10.1101/2022.09.15.508131>
4. <https://dx.doi.org/10.17504/protocols.io.5qpvo1nbl4o/v4>
5. <https://dx.doi.org/10.17504/protocols.io.eq2lyjzbrlx9/v2>
6. <https://community.artic.network/t/sars-cov-2-v4-1-update-for-omicron-variant/342>
7. <https://community.artic.network/t/sars-cov-2-version-5-3-2-scheme-release/462>

# Multi-pathogen & Syndromic Panels

“I think I know...”

---



# Hybrid Capture Enrichment Panels



01

Hybrid capture enables detection of a large range of potential pathogens and possible co-infections

02

More cost effective and sensitive than non-enriched shotgun metagenomic approaches<sup>1</sup>

03

Robust method can tolerate sequence differences and mismatches, allowing for successful enrichment of divergent sequences and variants.<sup>2</sup>

Panel	What it does
Viral Surveillance Panel (VSP)	Full genome of >60 viruses of high public health concern
Pan Coronavirus Panel (Pan CoV)	Full genome >200 coronaviruses and variants
Respiratory Virus Enrichment Kit (RVEK) (formerly RVOP)	Full genome of >40 common respiratory viruses
Respiratory Pathogen ID/AMR Panel (RPIP)	>280 respiratory pathogens (bacteria, viruses, fungi); >1200 AMR alleles; whole genome of SARS-CoV-2 and influenza
Urinary Pathogen ID/AMR Panel (UPIP)	>170 urogenital pathogens (bacteria, viruses, fungi, parasites) including STIs; >3700 AMR markers

<sup>1</sup> Bierstedt et al. Broad genomic AMR surveillance with an efficient, probe-captured NGS Panel. 2023 APHL poster

<sup>2</sup> Nagy-Szakal et al. Targeted hybridization capture of SARS-CoV-2 and metagenomics enables genetic variant discovery and nasal microbiome insights. <https://journals.asm.org/doi/10.1128/spectrum.00197-21>

# Importance of wastewater surveillance for public health

Community monitoring enables early detection of pathogens, threats and resistance profiles




Discover emerging microbes in entire populations irrespective of access to care



Detect & track pathogens, variants, and antimicrobial resistance (AMR) genes in communities early

## Urban monitoring of antimicrobial resistance during a COVID-19 surge through wastewater surveillance



Karthikeyan et al. 2021. Nature, <https://doi.org/10.1038/s41586-022-05049-6>

Harrington et al. 2022. Science of The Total Environment, <https://doi.org/10.1016/j.scitotenv.2022.158577>

# Comprehensive NGS approach for surveillance of viral and bacterial pathogens in wastewater

## Research Insights

Illumina's hybrid capture based Viral Surveillance Panel (VSP) and Urinary Pathogen ID/AMR Panel (UPIP) were evaluated on wastewater samples against traditional shotgun metagenomics


In the two extraction methods used, there was negligible difference in coverage and # of viral genomes

VSP enrichment of wastewater samples enhances (2x) the number of DNA and RNA viruses observed and allows for more granular strain identification

UPIP enrichment enhances detection, coverage and depth of urinary pathogens and AMR markers

AMR gene family coverage in UPIP enriched samples is 2.19X higher at 20X fewer reads

Broad surveillance of viral pathogens and AMR markers in wastewater, a complex environmental sample, at low sequencing depth



Janssen K, et al. Pathogen and antimicrobial resistance detection in wastewater using hybrid capture enrichment provides a method for tracking emerging infectious diseases and antimicrobial resistance genes. Poster presented at: Association of Public Health Laboratories Annual Conference; May 22-25, 2023; Sacramento, CA.

# Illumina Viral Surveillance Panel (VSP)

Culture-free panel for the detection and characterization of the most critical viral public health threats.



For Research Use Only. Not for use in diagnostic procedures.

## Applications

- Genomic Surveillance
- Seasonal Virus monitoring
- Emerging Virus Surveillance
- Wastewater Surveillance

## Methods

Probe-enhanced DNA and RNA Sequencing

## Other technologies

Sanger sequencing, qPCR (presence/absence), ELISA, pulse field electrophoresis (PFGE), Cell culture

## Features

Culture-free method

Broad viral surveillance panel including full genome design of 66 most critical viruses in Public Health, including, but not limited to; Mpox, SARS-CoV-2, Poliovirus, Influenza A/B, RSV, human parainfluenza viruses, and human adenoviruses.

Probe-enhanced panel, combined with Illumina RNA Prep with Enrichment Kit, generating reliable results

Analysis powered by DRAGEN Targeted Microbial App

Compatible with all Illumina platforms

## Current users:

- Researchers, Academic labs
- Public health labs
- Service providers

# Illumina Viral Surveillance Panel



## Arthropod transmitted / febrile tropical

Arbovirus  
**Chikungunya**  
**Dengue Virus 1**  
Dengue Virus 2  
Dengue Virus 3  
Dengue Virus 4  
**West Nile Virus**  
**Yellow fever virus**  
**Zika**



## Bloodborne

**Hepatitis A**  
Hepatitis B  
Hepatitis C  
Hepatitis E  
**HIV 1**  
HIV 2  
Torque Teno virus



## Encephalitis

Eastern equine encephalitis virus  
Hendra henipavirus  
Japanese encephalitis virus  
Nipah virus  
Tick-borne encephalitis virus  
Venezuelan equine encephalitis virus  
Western equine encephalitis virus



## Respiratory / Cardiopulmonary

Adenovirus  
Coronavirus -229E  
Coronavirus -HKU1  
Coronavirus- OC43  
Coronavirus-NL63  
**Influenza A**  
Influenza B  
Hantavirus  
MERS-CoV  
Metapneumovirus  
Parainfluenza  
Parechovirus  
**Respiratory syncytial virus**  
Rhinovirus  
SARS-COV  
**SARS-COV-2**



## Hemorrhagic fever

Chapare virus  
Crimean-congo hemorrhagic fever virus  
**Ebola**  
Guanarito Virus  
Rift Valley Fever Virus  
Junin Virus  
Kyasanur Forest disease virus  
**Lassa fever virus**  
Lujero hemorrhagic fever virus  
Machupo Virus  
**Marburg Virus**  
Omsk hemorrhagic fever virus  
Sabia virus



## Enteric

Aichivirus  
Astrovirus  
Coxsackievirus  
**Enterovirus**  
**Norovirus**  
**Poliovirus**  
Rotavirus  
Salivirus  
Sapovirus



## Oncolytic

**Oncolytic HPV**  
Polyomavirus



## Rash/Lesion

**Mpox**  
Parvovirus  
**Rubella**  
Variola virus

\*Sample types and sensitivity not validated by ILMN



# Pan-Coronavirus Panel (Pan-CoV)

Culture-free panel for the detection of known and closely related novel coronaviruses, for surveillance of animal reservoirs and human coronaviruses.



For Research Use Only. Not for use in diagnostic procedures.

## Applications

- Coronavirus surveillance
- Virus evolution & monitoring
- Coronavirus Identification
- Zoonotic, Spillover viral monitoring
- Wastewater surveillance

## Methods

Probe-enhanced (Hybrid-capture) Viral RNA

## Other technologies

Classic Microbiology Techniques, Sanger sequencing, qPCR (presence/absence), pulse field electrophoresis (PFGE)

## Features

Culture-free method

Whole-genome Sequencing of 203 known coronaviruses and over 370 strains of animal and closely related novel coronaviruses.

Enrichment panel allow for variation in viral genome, allowing for mutations and novel strains to be detected

Probe-enhanced panel for combining with Illumina RNA Prep with Enrichment Kit

Analysis powered by DRAGEN Targeted Microbial App

Compatible with all Illumina platforms



## Current users:

- Researchers, Academic labs
- Public health and Zoonotic labs
- Service providers

# Pan-Coronavirus Panel (Pan-CoV)

Monitor the environment and animal reservoirs to detect 'spill-over' to humans

## Illumina Pan-Coronavirus Panel (>200 coronaviruses)



Cui et al. Nature Reviews Microbiology 2019

# Illumina Respiratory Virus Enrichment Kit (RVEK)

A culture-free comprehensive respiratory viral panel to detect and characterize common respiratory viruses



For Research Use Only. Not for use in diagnostic procedures.

## Applications

- Genomic surveillance
- Seasonal Virus monitoring
- Respiratory Virus Identification
- Wastewater surveillance

## Methods

Probe-enhanced metagenomics (Precision Metagenomics)

## Other technologies

Classic virology techniques, Sanger sequencing, qPCR (presence/absence)

## Features

Culture-free method

Broad viral pathogen panel, including full genome design of >40 viruses

Robust method that remains effective with mutational changes that can render amplicon approaches or PCR ineffective

Includes human probes for positive internal control

Can detect as low as 2 genome copies per reaction

Analysis powered by DRAGEN Targeted Microbial App

Compatible with all Illumina platforms

## Current users:

- Researchers, Academic labs
- Applied Microbiology Labs
- (Pharma, Biotech)
- Service providers.

# Illumina Respiratory Virus Enrichment Kit



## Human adenovirus

Human adenovirus B1  
Human adenovirus C2  
Human adenovirus E4

## Human bocavirus

Human bocavirus 1 (Primate  
bocaparvovirus 1 isolate st2)  
Human bocavirus 2c PK isolate  
PK-5510  
Human bocavirus 3  
Human bocavirus 4 NI strain  
HBoV4-NI-385

## Human coronavirus

Human coronavirus 229E  
Human coronavirus HKU1  
Human coronavirus NL63  
Human coronavirus OC43  
Human coronavirus SARS-CoV-2

## Human enterovirus

Human enterovirus C104  
strain: AK11  
Human enterovirus C109  
isolate NICA08-4327

## Human metapneumovirus

Human metapneumovirus  
(CAN97-83)

## Human parainfluenza virus

Human parainfluenza virus 1  
Human parainfluenza virus 2  
Human parainfluenza virus 3  
Human parainfluenza virus 4a

## Human parechovirus

Human parechovirus 6  
Human parechovirus type 1  
PicoBank/HPev1/a

## Human respiratory syncytial virus

Human Respiratory syncytial virus  
(type A)  
Human Respiratory syncytial virus  
9320 (type B)

## Human rhinovirus

Human rhinovirus A89  
Human rhinovirus B14  
Human rhinovirus C (strain 024)

## Influenza A

Influenza A virus (A/Hong  
Kong/1073/99(H9N2))  
Influenza A virus  
(A/Korea/426/1968(H2N2))  
Influenza A virus  
(A/Michigan/45/2015(H1N1))  
Influenza A virus (A/New York/  
392/2004(H3N2))  
Influenza A virus (A/Puerto Rico/  
8/1934(H1N1))  
Influenza A virus  
(A/Texas/50/2012(H3N2))  
Influenza A virus (A/Zhejiang/  
DTID-ZJU01/2013(H7N9))  
Influenza A virus  
(AA/goose/Guangdong/1/1996  
(H5N1))

## Influenza B

Influenza B virus  
(B/Brisbane/60/2008)  
Influenza B virus  
(B/Colorado/06/2017)  
Influenza B virus (B/Lee/1940)  
Influenza B virus  
(B/Wisconsin/01/2010)  
Influenza B virus  
(B/Washington/02/2019)

## Polyomavirus

KI polyomavirus Stockholm 60  
WU Polyomavirus

## Human control genes

Human control genes

# Respiratory Pathogen ID/AMR Enrichment Panel Kit (RPIP)

Culture-free panel for identification, characterization, and quantitation of organisms associated with lower respiratory tract infections and genomic surveillance



For Research Use Only. Not for use in diagnostic procedures.

## Applications

- Antimicrobial Resistance
- Seasonal Virus monitoring
- Respiratory Pathogen Identification
- Genomic surveillance
- Wastewater surveillance

## Methods

hybrid capture enrichment  
(Precision Metagenomics)

## Other technologies

Classic Microbiology Techniques, Sanger sequencing, qPCR (presence/absence)

## Features

Culture-free method

Identify more than 282 pathogens associated with Lower Respiratory Tract Infection (52 Fungi, 42 Viruses, 182 Bacteria)

2,097 Antimicrobial Resistance Markers associated to 26 Antimicrobial Classes

Designed for the Identification of coinfections or hard-to-grow pathogens

Flexibility to accommodate a variety of sample types (e.g. sputum, bronchoalveolar lavage, wastewater, etc.)

Analysis powered by Explify RPIP Data Analysis

Compatible with all Illumina platforms

## Current user:

- Researchers, Academic labs
- Public health labs
- Clinical research
- Service providers

# Urinary Pathogen ID/AMR Enrichment Panel Kit (UPIP)

Culture-free identification and quantification of both common and underrecognized uropathogens that can cause urinary tract infections and sexually transmitted Infections.



For Research Use Only. Not for use in diagnostic procedures.

## Applications

- Antimicrobial Resistance
- Urinary Pathogen Identification
- Genomic surveillance
- Wastewater surveillance

## Methods

hybrid capture enrichment  
(Precision Metagenomics)

## Other technologies

Microbiology culture, Sanger sequencing, qPCR (presence/absence), STI Ag/Ab rapid test, shotgun metagenomics

## Features

Culture-free method

Identify more than 174 pathogens associated with Urinary Tract Infection (14 Fungi, 4 Parasites, 35 Viruses, 121 Bacteria)

Comprehensive Antimicrobial Resistance characterization of 3,761 Antimicrobial Resistance Markers associated to 18 Antimicrobial Classes

Flexibility to accommodate variations in sample types ( e.g., urine, swab, isolated organisms, cell culture and wastewater samples) and applications

Optimized library prep performance, generating reliable results

Automated analysis powered by Explify UPIP Data Analysis

Compatible with all Illumina platforms

## Current users:

- Researchers, Academic labs
- Public health labs
- Clinical research
- Service providers

# Illumina Infectious Disease Hybrid Capture Panel Comparison

	Illumina Viral Surveillance Panel VSP	Illumina Respiratory Oligo Virus Panel RVOP	Illumina Pan-Coronavirus Panel PCP	Illumina Respiratory Pathogen ID/AMR Panel RPIP	Illumina Urinary Pathogen ID/AMR Panel UPIP
<b>Target(s)</b>	66 viral pathogens, 203 strains	40 common respiratory viruses	225 Coronaviruses, over 370 strains	282 Respiratory pathogens (viral, bacteria, fungi) and 2108 antimicrobial resistance markers	174 Urinary pathogens (virus, bacteria, fungi, parasite) and 3728 antimicrobial resistance markers
<b>Whole Genome vs. Detection</b>	Whole Genome	Whole genome of common respiratory viruses, including recent flu strains and SARS-CoV-2. Includes human probes as a quality feature, human probes are included in every sample	Whole Genome	Detection of pathogens and AMR markers Whole genome: SARS-CoV-2, Influenza A/B	Detection of pathogens and AMR markers
<b>Library Prep</b>	RNA Prep with Enrichment	RNA Prep with Enrichment	RNA Prep with Enrichment	RNA Prep with Enrichment	DNA Prep with Enrichment
<b>Software</b>	DRAGEN Targeted Microbial App	DRAGEN Targeted Microbial App	DRAGEN Targeted Microbial App	Explify RPIP	Explify UPIP
<b>Configuration</b>	Standalone panel (96 reactions) Panel, Library Prep and Indexes Set A or B	Standalone panel	Standalone panel (96 reactions)	Panel, Library Prep, Indexes Set A, B, C, or D	Panel, Library Prep, Indexes Set A, B, C, or D



# Discovery

Hypothesis-free

---



# Pathogen discovery and whole-genome sequencing enables vaccine development driven by synthetic vaccinology

## SARS-CoV-2 Full Discovery

Shotgun sequencing for the discovery of the pathogen  
Clinical assay designed created once genome published

## SARS-CoV-2 Full Genome Sequence

The sequence was publicly available 31 days after the first case  
Vaccine development started immediately without isolation of viral particles

## Epitope & Vaccine Candidates

Bioinformatics screening for epitope candidates  
Package into vectors for animal test screening  
Worldwide Pharma, startups, and universities joined the race

## Clinical Trials

Moderna mRNA-based vaccine mRNA-1273 started  
clinical testing just 2 months after sequence identification

DISCOVERY TO WGS AND ANALYSIS IS THE FIRST STEP TO DISCOVER VACCINE ANTIGEN

# Illumina DNA Prep

Illumina library prep with the flexibility for whole-genome sequencing, amplicon & shotgun metagenomics applications for nucleotide resolution.



For Research Use Only. Not for use in diagnostic procedures.

## Applications

Genomic surveillance (HAI, Foodborne pathogens), gene-specific amplicons (16S, AMR genes), pathogen discovery, microbiome analysis & host/pathogen studies and more

## Methods

Whole genome sequencing, amplicon & shotgun metagenomics or metatranscriptomics

## Other technologies

Whole Transcriptome, qPCR.

## Features

Fastest Illumina library prep workflow, with ~3.5 hours total time

Flexibility to accommodate variations in sample types (gDNA, blood, saliva, dried blood spots, and direct from culture), DNA input amount (1-500 ng), and methods (WGS, amplicon or metagenomics)

Optimized library prep performance; delivers consistent insert sizes, uniform coverage, and optimized performance, regardless of DNA input amount or genome size.

The bead-based technology minimizes bias and opportunities for error, resulting in highly reproducible sequencing data.

Compatible with all Illumina platforms

Automatable through our many partners

## Current users:

- Researchers, Academic labs
- Public health labs
- Clinical research
- Applied Microbiology Labs (Pharma, manufacturing, direct to consumer)
- Service providers

# Ribo-Zero Plus Microbiome rRNA Depletion Kit

Fast, efficient depletion of undesirable host and pan-bacterial rRNA from complex microbial samples for metatranscriptomics research.



For Research Use Only. Not for use in diagnostic procedures.

## Applications

- Gut Metatranscriptomics
- Metatranscriptomics

## Methods

Metatranscriptomics,  
rRNA depletion

## Other technologies

Whole Transcriptome, qPCR.

## Features

Eliminate unwanted rRNA from bacteria common complex microbial samples, validated in adult and infant stool samples.

Optimized library prep performance, generating reliable results

Includes access to simplified data analysis and visualization with BaseSpace Microbiome Metatranscriptomics App

Compatible with MiSeq, NextSeq 550, NextSeq 1000, NextSeq 2000 and NovaSeq 6000.

## Current users

- Researchers, Academic labs
- Applied Microbiology Labs
- (Pharma, Biotechs)
- Service providers

# Infectious Diseases by the Numbers

770M	Cumulative cases of <b>COVID-19</b> <sup>1</sup>
400K	Respiratory deaths each year caused by <b>influenza virus</b> <sup>2</sup>
10.6M	People that fell ill from <b>tuberculosis</b> in 2021 <sup>3</sup>
39M	People living with <b>HIV</b> in 2022 <sup>4</sup>
3.8%	Global GDP lost due to <b>antimicrobial resistance</b> by 2050 <sup>5</sup>
247M	<b>Malaria</b> cases in 84 malaria-endemic countries in 2021 <sup>3</sup>
1.5M	New <b>hepatitis B</b> infections globally in 2019 <sup>3</sup>
\$28-45B	Direct annual cost of treating <b>healthcare associated infections</b> in the US <sup>6</sup>

Illumina Solution – building blocks to impact human health

COVIDSeq or IMAP

RVEK / RPIP or IMAP-Flu

Illumina and GenoScreen  
Deeplex Myc-TB Combo Kit

VSP or IMAP

RPIP and UPIP

IMAP

VSP

IDP

1. <https://covid19.who.int/> (Accessed Sept. 7, 2023); 2. <https://ourworldindata.org/influenza>; 3. World health statistics 2023: monitoring health for the SDGs, Sustainable Development Goals. Geneva: World Health Organization; 2023.; 4. <https://www.unaids.org/en/resources/fact-sheet>; 5. Global antimicrobial resistance and use surveillance system (GLASS) report 2022. Geneva: World Health Organization; 2022.; 6. Gidey K, et al. Clinical and economic burden of healthcare-associated infections: A prospective cohort study. PLoS One. 2023 Feb 23;18(2):e0282141.

# Thank you

**Tung, Vu**

09 June 2024| [tvu1@illumina.com](mailto:tvu1@illumina.com)